


INSTRUCTIONS FOR USE
ZILBRYSQ (ZIL-brisk)
(zilucoplan)
injection, for subcutaneous use
Single-Dose Prefilled Autoinjector



This Instructions for Use contains information on how to inject ZILBRYSQ.

Understanding Your ZILBRYSQ Prefilled Autoinjector

Read this Instructions for Use before using ZILBRYSQ autoinjector and each time you get a refill. There may be new information. This information does not take the place of talking to your healthcare provider about your medical condition or treatment. See the Medication Guide that comes with ZILBRYSQ for important information.

Your healthcare provider should show you how to prepare and inject ZILBRYSQ autoinjector properly before you use it for the first time. Do not inject yourself or someone else until you have been shown how to inject ZILBRYSQ correctly.

Keep this Instructions for Use and read it as needed until you have used all of the ZILBRYSQ autoinjectors in the packaging.

For general questions or help, please call your healthcare provider or UCBCares at 1-844-599-CARE (2273).

How Should I Store ZILBRYSQ Prefilled Autoinjectors?

- Store ZILBRYSQ autoinjectors in a refrigerator between 36°F to 46°F (2°C to 8°C) until the expiration date on the carton. **Do not freeze ZILBRYSQ.**
- ZILBRYSQ autoinjectors may be stored at room temperature up to 86°F (30°C) in the original carton for a single period of up to 3 months. Write the date you removed the autoinjectors from the refrigerator in the space provided on the carton. After ZILBRYSQ autoinjectors have been stored at room temperature, do not place them back in the refrigerator. Throw them away if not used within 3 months or if the expiration date has passed, whichever occurs first.
- Keep ZILBRYSQ autoinjectors in the original carton to protect them from light until it is time to use them.
- **Keep ZILBRYSQ autoinjectors and all medicines out of the reach of children.**

Important Information You Need to Know Before Injecting ZILBRYSQ:

- **Do not** accept or use ZILBRYSQ if the carton is damaged or opened or if the carton seals have been broken.
- **Do not** use ZILBRYSQ if the expiration date on the packaging has passed.
- **Do not** inject ZILBRYSQ more than 1 time per day.
- **Do not** use the ZILBRYSQ autoinjector if it has been dropped or looks damaged.
- **Do not** remove the needle cap from the ZILBRYSQ autoinjector until you are ready to inject.
- **Do not** use if the liquid has ever been frozen (even if thawed).
- **Do not** share or reuse the ZILBRYSQ autoinjector. The autoinjector is for 1-time (single-use) only. You may give or get an infection.
- **Do not** miss any doses of ZILBRYSQ. If you miss your ZILBRYSQ dose, inject a dose as soon as possible. Then, inject the next dose at your scheduled time. Do not inject more than 1 dose each day.

ZILBRYSQ Autoinjector Parts (Figure A):




Figure A

Preparing for Injection

Step 1: Take out the ZILBRYSQ prefilled autoinjector

If the ZILBRYSQ autoinjectors are stored in the refrigerator:

- Take the ZILBRYSQ autoinjectors carton out of the refrigerator and carefully remove 1 prefilled autoinjector from the carton. Put the rest of the autoinjectors in the carton back in the refrigerator.
- Before you inject ZILBRYSQ, let the autoinjector warm up to room temperature on a clean flat surface for **30 to 45 minutes**. This will help to reduce discomfort when injecting.
- **Do not** warm the ZILBRYSQ autoinjector in any other way (for example in a microwave, in hot water, or in direct sunlight).
- **Do not** remove the needle cap from the autoinjector until you are ready to inject.

If the ZILBRYSQ prefilled autoinjectors are stored at room temperature:

- Carefully remove 1 ZILBRYSQ prefilled autoinjector from the carton. Any remaining autoinjectors in the carton should not be placed in the refrigerator after it has been stored at room temperature.

Step 2: Inspect the ZILBRYSQ prefilled autoinjector

- Check the ZILBRYSQ autoinjector for damage (**Figure B**).
- **Do not** use if any part of the ZILBRYSQ autoinjector appears to be cracked, leaking, or broken.

- Check that the needle cap is not cracked or broken and that it's attached to the autoinjector (**Figure B**).
- **Do not** use if the needle cap is missing or is not securely attached.
- Check the expiration date and medicine name (ZILBRYSQ) on the autoinjector label (**Figure B**).
 - **Do not** use if the expiration date printed on the autoinjector has passed.
 - **Do not** use if the word ZILBRYSQ does not appear on the autoinjector label.
- Check the dose appearing on the label (**Figure B**) and make sure it matches your prescribed dose.
 - **Do not** use if the dose does not match your prescribed dose.
- Check the medicine through the Viewing window (**Figure B**). The medicine should be clear to almost clear and colorless. It is normal to see air bubbles.
 - **Do not** use if the medicine is cloudy, discolored, or has floating particles.




Figure B

Step 3: Gather supplies

Wash your hands with soap and water and dry them with a clean towel. Gather the following supplies on a clean, flat surface (**Figure C**):

- 1 ZILBRYSQ prefilled autoinjector
- 1 alcohol wipe (not supplied)
- 1 clean cotton ball or gauze pad (not supplied)
- 1 adhesive bandage (not supplied)
- 1 sharps disposal or puncture-resistant container (not supplied). See step 10 for instructions on throwing away the autoinjector.




Figure C

Step 4: Choose your injection site

Choose an injection site from the following areas:

- The stomach (abdomen), except for the 2-inch area around the belly button (navel) (**Figure D**).
- The front of the thighs (**Figure D**).




Figure D – For self or caregiver: abdomen and thighs

- The back of your upper arm (only if someone else is giving you the injection) (**Figure E**).




Figure E – For caregivers only: Back of upper arms

Choose a different injection site each time you give yourself an injection. If you want to use the same area of the body to inject, make sure it is at least 1-inch away from a spot you used before.

- **Do not** inject ZILBRYSQ into an area that is red, tender, bruised, swollen, hard or that has scars or stretch marks.

Step 5: Clean your injection site

- Clean the injection site with an alcohol wipe (**Figure F**).
- Let the skin dry for 10 seconds before injecting.
- **Do not** touch the injection site again before giving your injection.




Figure F

Injecting the Medicine

Step 6: Remove the needle cap

Hold the ZILBRYSQ autoinjector securely with one hand around the handle. Pull the needle cap straight off (away from you) with your other hand (**Figure G**).




Figure G

Throw away (discard) the needle cap into your household trash or a sharps disposal container.

Although you cannot see the needle tip, it is now uncovered.

- **Do not** touch the needle or let it touch anything.
- **Do not** touch the needle guard or put the needle cap back on. This could lead to a needle stick injury and could cause the autoinjector to release the medicine too soon or damage the needle inside.

Step 7: Hold the prefilled autoinjector

Hold the ZILBRYSQ autoinjector straight at a 90-degree angle to the cleaned injection site. Make sure that you can see the viewing window. **(Figure H).**




Figure H

Step 8: Inject the medicine

Place the ZILBRYSQ autoinjector straight against your skin, then firmly press the autoinjector down against your skin. **You will hear a first click.** This means the injection has started **(Figure I).**

- **Do not** lift the ZILBRYSQ autoinjector away from the skin during injection.
- **Keep holding the ZILBRYSQ autoinjector in place and pressed firmly against your skin until the viewing window is completely filled with the yellow indicator, and until you hear a second click.**
- The second click means that all the medicine has been injected and your ZILBRYSQ injection is now finished. You should see the yellow color indicator filling the Viewing window **(Figure J).**
- **Do not** remove the autoinjector away from the skin until you hear the second click because the medicine may not be fully injected.




Figure I




Figure J

Step 9: Remove the ZILBRYSQ autoinjector

- Remove the ZILBRYSQ autoinjector by carefully lifting the autoinjector straight up from your skin. The needle guard will automatically cover the needle. **Do not** touch the needle. **Do not** rub the injection site.
- **Important:** When you remove the autoinjector, if the Viewing window has not turned yellow this means you may not have received a full dose. Call your healthcare provider right away.

Step 10: Throw away (dispose of) the used ZILBRYSQ autoinjector

Throw away (dispose of) the used ZILBRYSQ autoinjector into an sharps disposal container (**Figure K**) right away. **Do not** throw away (dispose of) the ZILBRYSQ autoinjector in your household trash.

If you do not have an FDA-cleared sharps disposal container, you may use a household container that is:

- made of a heavy-duty plastic,
- can be closed with a tight-fitting, puncture-resistant lid, without sharps being able to come out,
- upright and stable during use,
- leak-resistant, and
- properly labelled to warn of hazardous waste inside the container.

Do not throw away (dispose of) your used sharps disposal container in your household trash unless your community guidelines permit this.

Do not recycle your used sharps disposal container.

Always keep the sharps disposal container out of the reach of children.




Figure K

Note: When your sharps disposal container is almost full, you will need to follow your community guidelines for the right way to dispose of your sharps disposal container. There may be state or local laws about how you should throw away your used needles and syringes. For more information about safe sharps disposal, and for specific information about sharps disposal in the state you live in, go to the FDA's website at:

<http://www.fda.gov/safesharpsdisposal>.

Step 11: Examine the injection site

Press a clean cotton ball or gauze pad over the injection site and hold it for 10 seconds (**Figure L**).

Do not rub the injection site. You may have slight bleeding. This is normal. Apply an adhesive bandage if needed (**Figure M**).




Figure L




Figure M

Your injection is now complete.

For more information or help, contact your healthcare provider.

If you have any questions, call UCBCares at 1-844-599-CARE (2273) or visit www.AskUCBCares.com.

Manufactured for:
UCB, Inc.
Smyrna, GA 30080

This Instructions for Use has been approved by the U.S. Food and Drug Administration

Issued: 10/2026